backstretcher®

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ENG: Read the instructions carefully before using the Neckstretcher bone.

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ENG: Lie on your back and bend your knee with your feet flat on the floor to release tension in your lower back.

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ENG: Place the neck stretcher under your head in such a way that the shoulders rest on the long part with the three rollers where contact is made with the tense muscles. The upper roller is in contact with the pressure points at the base of your skull and both sides of your cervical vertebrae.

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ENG: After taking this position, ease yourself by slowly lowering your legs.

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ENG: Always lie down slowly and quietly and never force the movement.

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ENG: Start sessions that last no more than a few minutes. As you progress, gradually build this up minute by minute for a maximum of 10-15 min.

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Important: If you have vertebrae that have grown together (so-called block vertebrae), it is wise to consult your doctor or therapist before using the neck stretcher. We also recommend people with severe osteoporosis, pregnant women and those who have other doubts about suitability to discuss with their doctor or therapist whether the use of the neck stretcher can also be beneficial in their case.

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